

Matcher Version 4.5 – April 2017

Overview – what is this release for?

This release provides enhancements to functionality in a number of areas. We have a new cryo label sheet available that includes PCR tube labels for PGS/PGD. There is a new customisable patient list as well as changes around the way you open new cycles to make the system more intuitive. We have made some changes to the No Match result screen to improve auditability and make 'no matches' clearer. We have also implemented updated functionality to allow partners to be imported from a linked Patient Management System where the patient is already in Matcher (and vice versa). There are also several other changes and minor fixes which are detailed below.


This is an update to the Matcher PC application, Pocket Matcher application and database.

What is included in this release?

1. [PGS labelling](#)
2. [Updated patient list](#)
3. [Updated cycle types](#)
4. [Updated 'No Match' screen](#)
5. [Enhancements to patient import facility for patient couples](#)
6. [Updated Pocket Matcher witnessing images](#)
7. [Custom patient ID field on Pocket Matcher](#)
8. [Improved error handling](#)
9. [Ability to archive / remove 'Matcher Types' \(drop-down list options\)](#)
10. [New override reports](#)
11. [Various other small changes and minor fixes](#)

1. PGS labelling


Developed with PGS labs, we have a new cryo label sheet available containing PCR tube labels (in addition to the standard straw and vial labels). These curved labels (to fit tapered PCR tubes) are suitable for 0.2µl and 0.5µl tubes and are self-laminating. The labels contain a 2d barcode and patient / item information. These labels are already in use in PGS labs to identify samples throughout the PGS process.



2. Updated patient list

We have made some changes to the patient list screen to allow you to **customise the information you see** – you can view information such as the **last procedure performed, the next procedure(s) due, the number of days since the cycle started (i.e. number of days since the first procedure)** etc. You can also remove information that you don't need.

We have set the default view to show the columns that we think would be most useful for you (see below), but you can add or remove columns as you wish.



Matcher

File Reports Tools Help

User: Matcher Admin 07

Patient Search

Enter Name or Clinic Id Clear Search

☒ All Patients ☐ Active Cycles

Patient ID	National ID	Patient Name	Latest Cycle	Cycle Day	Last Procedure	Next Expected Procedure(s)	Process Map(s)	Cycle Name
23580	454-046-286	Mountie, Rocky	2017	4	ID to dish	Dish to dish	Oocyte recovery	IVF / ICSI (Major)
95938457-C	PP		7	4	Dish to dish	Dish to wash dish	Oocyte recovery	IVF / ICSI (Major)
25648	CS		7	3	Gradient tube to wash	Wash dish to dish	IVF	Insemination, Frozen Embryo Transfer (Minor)
TEST000003	123		7		Oocyte recovery	Dish to dish	IVF	IVF / ICSI (Major)
TEST000004	180379-SSSS	TEST MALE DONOR, Test Male						
TEST000002	987654321098	TEST PARTNER, Test Partner						
23574	ABC1234	Flynn, Gemma						
23575	ABC1235	McCaw, Richie						
23576	999999998	Sheridan, Andrew						
23582		Shepherdson, Guy						
23579	046-454-286	Mountie, Appalachia						

Log Off + Add → Select

Cycles remaining: 234 Major, 248 Minor, 250 Other | 0 5 | File Version: 4.5.6312.27559 | Database Version: 4.5.0.0

To add or remove columns, right click on the patient list as shown on the following page.

Patient ID	National ID	Patient Name
0124545		Clarke, Hannal
16543	123456789012	Hildebrandt, L
7897897895		
134652346		
201611221100		
555444667		
0123245984		
777888899		
201610111211		
2356487		
U-Y-BL-120585-2-F		

Reset Patient Sorting

Allow Change To Column Order

- ☒ National ID
- ☒ Patient Name
- ☐ Donor
- ☐ Partner Matching
- ☐ Donor IDs
- ☒ Latest Cycle
- ☒ Cycle Day
- ☐ Cycle Number
- ☒ Cycle Name
- ☐ Cycle Description
- ☐ Cycle Finished
- ☐ First Procedure Date
- ☒ Last Procedure
- ☐ Last Procedure Date
- ☒ Process Maps
- ☐ Process Map Routes
- ☐ Expected Procedure Types
- ☒ Next Expected Procedures
- ☐ Partner ID
- ☐ Partner Name

Right-click (or long press on a touchscreen) to bring up this menu.

Click 'Reset Patient Sorting' to order the list by Latest Cycle, Cycle Type, Patient Name.

Click 'Allow Change To Column Order' to allow the columns in the list to be moved.

All columns currently shown in the patient list are ticked.

Click on any ticked column names to remove them from the patient list.

Click on any un-ticked columns to add them to the patient list.

Any changes you make to the patient list are specific to each PC/Windows user account so you can have a different view in the office where you may need to see different information compared to the lab.

An alternative layout, which may be useful for a reception or administration machine showing whether the patient is a donor and also the patient's partner details is shown on the next page:

Matcher

File Reports Tools Help

User: Matcher Admin 09

Patient Search

Clear Search

☒ All Patients ☐ Active Cycles

Patient ID	National ID	Patient Name	Donor	Partner ID	Partner Name
0123456789		Adams, Amy	No		
U-Y-C-060790-1570-F	120456784612	Art, Adam	No		
201611221100	4561245878	Bailey, Sara	No		
555444667	1597897978	Blake, Emily	No	789452548	Ponting, Tom
ROB123	223222222	Bob, Rob	No		
0123245984		Bond, Emma	No	124585	Leto, Matthew
0124545		Clarke, Hannah	Yes		
12327894	262346246245	Clarke, Gabriella	No		
134652346	234634574567	Cycle, Test	Yes	123555444	Hussey, Chantelle
4567456	45874568	Cycler, Test1	No	201610111211	Paul, Steff
65465454		Finch, Rebecca	No		
45678945	789789789	France, Augusta	No		
16543	123456789012	Hildebrandt, Leonie	No	16544	Hildebrandt, Lukas
16544	987654321089	Hildebrandt, Lukas	No	16543	Hildebrandt, Leonie
123555444		Hussey, Chantelle	No	134652346	Cycle, Test
58578784		Jackson, Benedito	No		
15151559	852741963	Jerry, Ben	No		
201610141045		Lee, Rob	Yes		

Log Off

Cycles remaining: 9920 Major, 9985 Minor, 9984 Other

0 5

File Version: 4.5.6320.22749 Database Version: 4.5.0.0

As you can now see more 'live' data on this screen, we have set it to refresh every minute, but you can increase this time period if you wish (via **Tools > Settings > Patient List Refresh Rate**).

There is also a small update to the Pocket Matcher patient list whereby you will now see patients with the most recently added cycles at the top of the list.

3. Updated cycle types and descriptions

We have updated the cycle types available to make the system more intuitive and easy to use. The cycle types now match the terms that you use in the clinic, such as IVF/ICSI. The new cycle types available for selection are shown below:

Cycle Type Selection

- IVF / ICSI (Major)
- Frozen Embryo Transfer (Minor)
- IUI / Donor Insemination (Minor)
- Elective / Donor Sperm Cryopreservation (Minor)
- Elective / Donor Oocyte Cryopreservation (Minor)
- Pre-Cycle Semen Assessment (Other)
- Pre-Cycle Blood Testing (Other)

Description of selected cycle

IVF or ICSI, which may include IMSI, GIFT, gamete thawing, altruistic fresh oocyte sharing, fresh embryo transfer, and/or cryopreservation of excess material at the end of the cycle

OK Cancel

When you select a cycle type, you will also see a further description as shown on the left.



Good to know...

The pricing of cycles is not changing. This change is to make it easier to understand which type of cycle to select using terminology most familiar to you.


If we have previously set up different cycle types for you, we will ensure these are kept in place.

4. Updated 'No Match' result screen

We have made some changes to the 'No Match' result screen, both to make the 'No Match' result clearer and to improve auditability.



In the event of a No Match result, you will now be asked to select a reason for the No Match. This makes it easier for anyone using touchscreens in particular to select a reason from a pre-defined list, and helps you with root-cause analysis.

We have also updated the screen to make the No Match more obvious.



This list, as with the contents of all Matcher drop-down lists, is configurable to your clinic so that you can define your own reasons to fit in with the way that you work, although we provide an out-of-the-box list that covers the main reasons.

If you are set up as an administrator, you can update the list for your own clinic. To do this, go to **Tools > Matcher Type Management > No Match Reason** as shown on the next page.



5. Enhancements to import facility for patient couples

We have updated the import facility to allow you import a patient or partner from a source patient management system where one of the couple is already in the Matcher database. Previously, if one of the couple was already in Matcher, the other would have to be added manually.

6. Updated Pocket Matcher witnessing images

We noticed that in some cases the text stamp on Pocket Matcher witnessing images could potentially obscure some of the image. We have therefore updated the images to maximise the image space. This is shown below:

Old




New



7. Custom patient ID field on Pocket Matcher


Where you have the custom patient ID field enabled for patients (this is an additional field that some clinics use to hold a National ID, Passport number or similar), this is now available to view and search on the Pocket Matcher.

If you do not use this field and would like to, please see the Matcher User Guide, or contact our Support Team who will be happy to help you enable the field for use.





8. Improved error handling

In the unlikely event that you encounter an error when using Matcher, you will now see a friendly message, rather than a system message.



9. Ability to archive / remove 'Matcher Types' (drop-down list options)

The values in all drop-down lists in Matcher can be controlled through Matcher Types (accessible through **Tools > Matcher Type Management** for administrators) allowing you to customise Matcher for your own clinic. We have now added further flexibility by allowing these values to be archived, meaning the value can no longer be selected, but any previous selections of that value remain, and can be viewed in reports and summary lists. See diagram on the next page.



10. New override reports

We had become aware of a few scenarios in which the override reports (for a specific patient/cycle and all patients) were not displaying as expected. We have updated these reports with new versions displaying correctly in all scenarios and updated the layouts slightly to improve usability.

11. Various small changes and bug fixes

- (a) Fixed various minor label report layout issues.
- (b) Fixed various issues and enhanced functionality on the **Patient > Cryo** page. For example, you can now add multiple batches/lots of items.
- (c) Updated functionality to ensure all reports have descriptive tooltips.
- (d) Fix to Batch/Lot Summary Report to ensure images always displayed.
- (e) Enhanced functionality to ensure new batch/lot items default to correct category.
- (f) Standardised placement of buttons, including Close and Delete buttons.
- (g) Fixed an issue with the process map search facility.
- (h) Added specific Administrator users for all IMT support staff to ensure full auditability.
- (i) Fixed an issue with the camera export / import facility which will allow any camera issues to be resolved more easily.
- (j) Updated all label report names starting 015 to 017 to correlate to the updated cryo label sheets.
- (k) Consolidated the ID card reports to make implementation easier plus a small update to the layout to accommodate longer ID numbers.
- (l) Added the facility for a small delay in the Patient search facility. This can improve performance and speed of the search on slow networks. This is switched off by default (set to 0), but can be switched on through **Tools > Settings > Patient List Search Delay** and specifying the number of seconds' delay required. We would recommend setting this somewhere between 0 and 3 seconds.
- (m) Bug fix to ensure ID card images are always presented at the correct size.
- (n) Bug fix to ensure the patient custom field is not searchable if not enabled.
- (o) Bug fix to ensure you cannot save a process map without a type selected.
- (p) Updated the layout for Matcher Types to allow more text to be entered for the description.

What do I need to do to implement this update?

You do not need to do anything. Our Support Team will contact your IT Team to arrange installation / upgrade. This can now be done with our one-step installer making upgrades much quicker and easier.

What else do I need to do?

You do not need to do anything further.

You can contact the Matcher Support Team on +44 (0)1829 771 327 or by email at matcher@imtinternational.com if you have any questions.

Documentation updates

The User Guide has been updated with this release. You can find this by clicking **Help > User Guide** from within Matcher.

Are there any impacts or risks in applying this update?

- There should not be any disruption to your ongoing work during the upgrade, but access to each PC running the Matcher software will be required to perform the upgrade and user acceptance testing, during which it will be unavailable to use.
- This is a low risk change and does not make any changes to patient or witnessing data. The changes have undergone thorough system and acceptance testing at IMT
- In the unlikely event that the system upgrade causes any negative issues, the software can be rolled back to a previous version.

Quality assurance (QA)

These changes have been subject to IMT's quality assurance process as follows:

- All system changes are logged, tracked and source controlled.
- System testing is completed throughout the development lifecycle to ensure changes have been completed as required.

- Full regression testing is completed to ensure changes have had no adverse impacts on other areas of functionality.
- Full internal 'user acceptance testing' is completed by members of the IMT team.
- Usability testing is incorporated into all stages of the QA process.
- Beta testing is completed by a limited number of Matcher customers at their sites with live data (normally one to three customers, which are selected based on the suitability of their workflows with the areas of affected functionality).

Further questions

If you have any questions on this release, please contact the Matcher Support Team on +44 (0)1829 771 327 or email matcher@imtinternational.com.