

## Matcher Version 4.5.1 – May 2017 Release Notes

### Overview – what is this release for?


This release provides some minor updates to the patient list and the patient matches screens based on feedback received from customer testing of version 4.5.

This is an update to the Matcher PC application only.

### What is included in this release?

#### 1. Added 'No Match Reason' to patient matches screen

We have added the 'No Match Reason' to the patient matches screen as shown below:



ID	Match Date	User	Match Result	Procedure	Notes
625	28/04/2017 16:19	Matcher Admin 09	Wrong Patient - Wrong patient selected	OPU Tube to Egg Search Dish	
624	28/04/2017 16:18	Matcher Admin 09	Wrong Patient - Wrong dish / tube / item used	OPU Tube to Egg Search Dish	
623	28/04/2017 15:54	Matcher Admin 09	Matched	Female ID to OPU Tube	

Where you have a 'No Match' or 'Wrong Patient' result and have recorded a 'No Match Reason' it will show in the 'Match Result' column.

#### 2. Changes to patient list

We've made some minor changes to the patient list screen to add further functionality, as follows:

- You can right-click and choose **Allow Change To Column Size** which will allow you to change the column size. This means that you cannot accidentally change the column size without first selecting this option.
- You can right-click and choose **Reset To Column Defaults** if you want to change the columns back to the Matcher default columns, which are:
  - Patient ID
  - Patient Custom Field (if used)
  - Patient Name
  - Latest Cycle
  - Cycle Day
  - Last Procedure
  - Next Expected Procedure(s)
  - Active Process Map(s)
  - Cycle Name



Patient Name		Last Procedure	
Mountie, Rocky			
Hildebrandt		Reset Patient Sorting	
Sharapova	<input checked="" type="checkbox"/>	Reset To Column Defaults	
Norton, He	<input checked="" type="checkbox"/>	Allow Change To Column Order	
TEST PATI	<input checked="" type="checkbox"/>	Allow Change To Column Size	
Perez Rod	<input checked="" type="checkbox"/>	National ID	
María Rosa	<input checked="" type="checkbox"/>	Patient Name	
		Donor	
		Partner Matching	

There are two new options available when you right-click on the patient list:

- 'Reset To Column Defaults' resets the columns back to the default.
- 'Allow Change To Column Size', when ticked allows you to change the size of the columns.

Note that it is possible to re-size columns on the patient list using a touchscreen monitor or tablet, but it is more difficult. If you are having issues resizing columns, you may want to temporarily use a mouse whilst you set up the list to your preference.

### What do I need to do to implement this update?

You do not need to do anything. Our Support Team will contact your IT Team to arrange installation / upgrade. This can now be done with our one-step installer making upgrades much quicker and easier.

### What else do I need to do?

You do not need to do anything further.

You can contact the Matcher Support Team on +44 (0)1829 771 327 or by email at [matcher@imtinternational.com](mailto:matcher@imtinternational.com) if you have any questions.

### Documentation updates

The User Guide has been updated with this release. You can find this by clicking **Help > User Guide** in Matcher.

### Are there any impacts or risks in applying this update?

- There should not be any disruption to your ongoing work during the upgrade, but access to each PC running the Matcher software will be required to perform the upgrade and user acceptance testing, during which it will be unavailable to use.
- This is a low risk change consisting of minor changes to the PC application. It does not affect the Matcher database, nor any data.
- In the unlikely event that the system upgrade causes any negative issues, the software can be rolled back to a previous version.

### Quality assurance (QA)

These changes have been subject to IMT's quality assurance process as follows:

- All system changes are logged, tracked and source controlled.
- System testing is completed throughout the development lifecycle to ensure changes have been completed as required.
- Full regression testing is completed to ensure changes have had no adverse impacts on other areas of functionality.
- Full internal 'user acceptance testing' is completed by members of the IMT team.
- Usability testing is incorporated into all stages of the QA process.
- Beta testing is completed by a limited number of Matcher customers at their sites with live data (normally one to three customers, which are selected based on the suitability of their workflows with the areas of affected functionality).

### Further questions

If you have any questions on this release, please contact the Matcher Support Team on +44 (0)1829 771 327 or email [matcher@imtinternational.com](mailto:matcher@imtinternational.com).